

<u>Transport of watershed-based materials: Sediment studies:</u> (A. Winter, N. Carballeira)

Assess historical rates of cross-shelf transport by determining relative quantities of terrestrial/marine materials in sediments. Sediment cores will be taken at three nearshore/offshore sites and processed by standard techniques (photography, grain size analysis, organic content, etc.). Terrestrial biomarkers (δ^{13} C, C₂₅-C₃₅ n-alkanes, and C₂₃- C_{34} fatty acids, etc.) and marine biomarkers (C_{17} - C_{20} n-alkanes, cholesterol, dinosterol, etc.) will be analyzed following techniques of Ohkouchi et al. (1997). Pb-210/Cs-137 and non-destructive dating techniques will also be used to determine age and mass accumulation rates.


Sections of the sediment samples analyzed (LP 5,6,12)





An increase of sand with depth

Increase in terrestrial sediment supply



A typical shelf ²¹⁰Pb profile

Fatty Acid Analysis of Sediments


- Biological origin of sedimentary organic matter possible biological mixing
- Assignment of specific fatty acids as markers for the contribution of particular biota

Isolation Procedure and Characterization


Total lipids extracted from sediment with CHCl₃/MeOH

- 1. GC-MS analysis of the total lipids extracted
- 2. Transesterification with HCI/MeOH
- 3. GC-MS analysis of the total fatty acid methyl esters


Fatty Acid Methyl Esters (LP12)




Fatty Acid Abundance as a Function of Depth (cm) for LP12




Pb-210 Activity (dpm/g)



Fatty Acid Abundance as a Function of Depth (cm) for LP5




Pb-210 Activity (dpm/g)




The profiles are relatively linear and likely provide a reasonable estimate on the modern rate of sediment accumulation; biological and physical mixing in this environment is probably <10 cm.

Fatty Acid Abundance as a Function of Depth (cm) LP6



Pb-210 Activity (dpm/g)



Fatty Acid Abundance as a Function of Depth (cm) LP3

